

DIMA PRINT STONE GRAY

THE COLOR FOR CERAMISTS

A NEUTRAL BACKDROP AND TOP AESTHETICS



MORE INFO
[KULZER.COM/
STONEGRAY](http://KULZER.COM/STONEGRAY)



KULZER
MITSUI CHEMICALS GROUP

DIMA PRINT STONE GRAY

THE AESTHETIC SPECIALIST

- ▶ **Traditional look:** Ensures natural color perception and no visual discoloration for tasks that need top aesthetics
- ▶ **Traditional feel:** Doesn't distract; feels like gray gypsum when handling, grinding, etc.
- ▶ **Traditional properties:** No changes to how you work with the model thanks to its "non-stick" surface and high heat resistance



Printing performance overview

Example	One lower jaw full (hollow)
Build height [mm]	27.9
Time [min] for 50 µm resolution	61
Time [min] for 70 µm resolution	41
Time [min] for 100 µm resolution	29
Material consumption [ml]	22.5
Average material costs [€]	6.16

Product overview

Color	gray
Cleaning time in Isopropanol (cara Print Clean or ultra sonic bath)	Pre-cleaning: 3 minutes Post-cleaning: 2 minutes
Post-curing time (HiLite Power 3D)	10 minutes (5 minutes front side + 5 minutes reverse side)
Post-curing time LEDcure (no turning of object necessary)	Select dima Print Stone gray program

Overview of dima Print Stone materials

Stone beige The all-rounder	Stone teal The ortho specialist	Stone gray The aesthetic specialist

Part of the validated 3D-Printing pro Solution workflow



Watch video!
[kulzer.com/3dp-video](https://www.kulzer.com/3dp-video)

